




Fast Detection · High Sensitivity



CISH analysis of a tonsil tissue using the ZytoFast[®] human Ig-kappa/Ig-lambda Permanent CISH Kit.




ZytoFast[®] EBV Probe hybridized on EBV infected tonsil tissue.




ZytoFast[®] HPV High-Risk Types Probe on condyloma tissue with HPV positive cells.


DAB



Permanent Red



NBT/BCIP



Selected Products of the ZytoFast[®] PLUS Portfolio

For more product information please contact info@zytovision.com or your local dealer.

Prod. No.	Product	Label	Tests* (Volume)
T-1113-400	ZytoFast CMV Probe	Digoxigenin	40 (400 µl)
T-1114-400	ZytoFast EBV Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1055-400	ZytoFast HPV type 6/11 Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1056-400	ZytoFast HPV type 16/18 Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1057-400	ZytoFast HPV type 31/33 Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1140-400	ZytoFast HPV High-Risk (HR) Types Probe (specific for HPV type 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82) C€ IVD	Digoxigenin	40 (400 µl)
T-1144-400	ZytoFast HPV Screening Probe (specific for HPV type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82) C€ IVD	Digoxigenin	40 (400 µl)
T-1017-400	ZytoFast human Ig-kappa/Ig-lambda Probe C€ IVD	Digoxigenin + Biotin	40 (400 µl)
T-1115-400	ZytoFast human Ig-kappa Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1116-400	ZytoFast human Ig-lambda Probe C€ IVD	Digoxigenin	40 (400 µl)

Related Products

Complete Kits

Prod. No.	Product	Label	Tests* (Volume)
T-1005-40	ZytoFast human Ig-kappa/Ig-lambda CISH Kit C€ IVD Incl. Ig-kappa/Ig-lambda Probe (Digoxigenin/Biotin labeled), 0.4 ml; 28S rRNA (+) Control Probe (Biotin labeled), 0.1 ml; RNA (-) Control Probe (Biotin labeled), 0.1 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-Biotin/DIG-Mix, 4 ml; AEC Solution, 4 ml; NBT/BCIP Solution, 4 ml; Nuclear Green Solution, 20 ml; Mounting Solution (aqueous), 4 ml	Digoxigenin + Biotin	40
T-1105-40	ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit C€ IVD Incl. Ig-kappa/Ig-lambda Probe (Digoxigenin/Biotin labeled), 0.4 ml; 28S rRNA (+) Control Probe (Biotin labeled), 0.1 ml; RNA (-) Control Probe (Biotin labeled), 0.1 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 2x 50 ml; Anti-Biotin/DIG-Mix, 4 ml; HRP-Green-Solution A, 0.8 ml; HRP-Green-Solution B, 15 ml; Permanent Red Solution A, 0.25 ml; Permanent Red Solution B, 15 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	Digoxigenin + Biotin	40

For the Detection of Digoxigenin-labeled ZytoFast[®] Probes

T-1061-40	ZytoFast PLUS CISH Implementation Kit AP-NBT/BCIP C€ IVD Incl. DNA (+) Control Probe, 0.1 ml; DNA (-) Control Probe, 0.1 ml; 28S rRNA (+) Control Probe, 0.1 ml; RNA (-) Control Probe, 0.1 ml; Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x 50 ml; Rabbit-anti-DIG, 4 ml; Anti-Rabbit-AP-Polymer, 4 ml; NBT/BCIP Solution, 4ml; Nuclear Red Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	Digoxigenin	40
T-1063-40	ZytoFast PLUS CISH Implementation Kit HRP-DAB C€ IVD Incl. DNA (+) Control Probe, 0.1 ml; DNA (-) Control Probe, 0.1 ml; 28S rRNA (+) Control Probe, 0.1 ml; RNA (-) Control Probe, 0.1 ml; Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x 50 ml; Mouse-anti-DIG, 4 ml; Anti-Mouse-HRP-Polymer, 4 ml; DAB Solution A, 0.3 ml; DAB Solution B, 10 ml; Nuclear Blue Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	Digoxigenin	40
T-1151-40	ZytoFast PLUS CISH Implementation Kit AP-Permanent Red C€ IVD Incl. DNA (+) Control Probe, 0.1 ml; DNA (-) Control Probe, 0.1 ml; 28S rRNA (+) Control Probe, 0.1 ml; RNA (-) Control Probe, 0.1 ml; Heat Pretreatment Solution EDTA, 500 ml; Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 4x 50 ml; Rabbit-anti-DIG, 4 ml; Anti-Rabbit-AP-Polymer, 4 ml; Permanent Red Solution A, 0.25 ml; Permanent Red Solution B, 15 ml; Mayer's Hematoxylin Solution, 20 ml; Mounting Solution (alcoholic), 4 ml	Digoxigenin	40

Digoxigenin-labeled Control Probes

T-1053-400	ZytoFast DNA (+) Control Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1054-400	ZytoFast DNA (-) Control Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1120-400	ZytoFast 28S rRNA (+) Control Probe C€ IVD	Digoxigenin	40 (400 µl)
T-1119-400	ZytoFast RNA (-) Control Probe C€ IVD	Digoxigenin	40 (400 µl)

* Using 10 µl probe solution per test. C€ IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.

