

Indication Overview

Solid Tumors

Brain and Neural Tumors
Breast Cancer
Cervical Cancer
Lung Cancer
Prostate Cancer
Renal Cell Carcinoma
Sarcoma

Hematology

Acute Lymphoblastic Leukemia (ALL)
Acute Myelogenous Leukemia (AML)
Chronic Lymphocytic Leukemia (CLL)
Chronic Myelogenous Leukemia (CML)
Multiple Myeloma
Myelodysplastic Syndrome (MDS)
Lymphoma

Genetics

Sex Mismatched Bone-Marrow Transplantant Management
Prenatal, Postnatal, and Preimplantation Genetics

Selected Products of the ZytoLight® Portfolio

For more product information please contact info@zytovision.com or your local dealer.

Indication	Prod. No.	Product	Label	Tests* (Volume)
Brain and Neural Tumors	Z-2075-50/-200	ZytoLight SPEC 1p36/1q25 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2076-50/-200	ZytoLight SPEC 19q13/19p13 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2063-200	ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe CE IVD	●/●	20 (200 µl)
Breast Cancer	Z-2198-50	ZytoLight SPEC TP53/17q22 Dual Color CE IVD	●/●	5 (50 µl)
	Z-2015-50/-200	ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2190-50	ZytoLight SPEC ERBB2/D17S122 Dual Color Probe CE IVD	●/●	5 (50 µl)
	Z-2069-50/-200	ZytoLight SPEC ESR1/CEN 6 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2072-50/-200	ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
Lung cancer	Z-2117-50/-200	ZytoLight SPEC ALK/EML4 TriCheck Probe CE IVD	●/●/●	5/20 (50/200 µl)
	Z-2124-50/-200	ZytoLight SPEC ALK Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2087-200	ZytoLight SPEC MET/CEN 7 Dual Color Probe CE IVD	●/●	20 (200 µl)
	Z-2194-200	ZytoLight SPEC NRG1/CD74 TriCheck Probe CE IVD	●/●/●	20 (200 µl)
	Z-2179-200	ZytoLight SPEC CD274, PDCD1LG2/CEN 9 Dual Color Probe CE IVD	●/●	20 (200 µl)
	Z-2148-50/-200	ZytoLight SPEC RET Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2144-50/-200	ZytoLight SPEC ROS1 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
Sarcoma	Z-2116-50/-200	ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2100-50/-200	ZytoLight SPEC DDT3 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2096-50/-200	ZytoLight SPEC EWSR1 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2183-50	ZytoLight SPEC EWSR1/FLI1 TriCheck Probe CE IVD	●/●/●	5 (50 µl)
	Z-2139-50	ZytoLight SPEC FOXO1 Dual Color Break Apart Probe CE IVD	●/●	5 (50 µl)
	Z-2185-50	ZytoLight SPEC FOXO1/PAX3 TriCheck Probe CE IVD	●/●/●	5 (50 µl)
	Z-2132-50	ZytoLight SPEC JAZF1 Dual Color Break Apart Probe CE IVD	●/●	5 (50 µl)
	Z-2013-50/-200	ZytoLight SPEC MDM2/CEN 12 Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2097-50/-200	ZytoLight SPEC SS18 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2184-50	ZytoLight SPEC SS18/SSX1 TriCheck Probe CE IVD	●/●/●	5 (50 µl)
	Z-2109-50/-200	ZytoLight SPEC TFE3 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
	Z-2151-50	ZytoLight SPEC USP6 Dual Color Break Apart Probe CE IVD	●/●	5 (50 µl)
	Hematology	Z-2216-50	ZytoLight SPEC 11q gain/loss Triple Color Probe CE IVD NEW	●/●/●
Z-2192-50/-200		ZytoLight SPEC BCL2 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
Z-2177-50/-200		ZytoLight SPEC BCL6 Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
Z-2111-50/-200		ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe CE IVD	●/●	5/20 (50/200 µl)
Z-2160-50/-200		ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe CE IVD	●/●/●	5/20 (50/200 µl)
Z-2110-50/-200		ZytoLight SPEC IGH Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
Z-2210-50		ZytoLight SPEC IRF4, DUSP22 Dual Color Break Apart Probe CE IVD NEW	●/●	5 (50 µl)
Z-2090-50/-200		ZytoLight SPEC MYC Dual Color Break Apart Probe CE IVD	●/●	5/20 (50/200 µl)
Z-2209-50		ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Probe CE IVD NEW	●/●/●	5 (50 µl)
Z-2159-50/-200		ZytoLight SPEC TP53/ATM Dual Color Probe CE IVD	●/●	5/20 (50/200 µl)

Related Products

Z-2028	ZytoLight FISH-Tissue Implementation Kit CE IVD	5/20
<small>5: Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml 20: Incl. Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 500 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml</small>		
Z-2099	ZytoLight FISH-Cytology Implementation Kit CE IVD	20
<small>Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl₂, 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTect-Solution, 0.8 ml</small>		

* Using 10 µl probe solution per test. CE IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.

